

Experience the beauty the perfect The Chart Projector.



EUCARIS **tSCP**

800
Chart Projector

The TSCP-800 Chart projector provides clear and sharp images with fast chart swapping.

Features of TSCP-800

When operate TSCP-800, You can expect a high quality and best customer service.

Simple, Clear and sharp image

The TSCP-800 Chart projector provides clear and sharp images with fast chart swapping.

High Quality

The TSCP-800 has high resolution of 50 lines at full fields.

Capability of various test

The TSCP-800 with Red/Green and polarizing filters are for binocular balance test, stereo test, aniseikonia test, fusion test.

Specification (TSCP-800)

Refraction Distance	2.5 ~ 8m
Projection Distance	1.5 ~ 10m
Projection Magnification	30x (in 5m refraction)
Chart Rotation Speed	1 frame Average 0.2 sec
Number of Chart	32
Power Supply	AC 100-240V 50/60Hz, 100VA
Lamp	LED 10W
Power Consumption	0.2A at 220V
Power Saving	3min
Tilt Angle	±20°
Dimension	310(D) X 170(W) 230(H) mm
Weight	3.9 kg
Accessories	Remote Control / Table Stand / Screen / Dust Cover / Fuse / Power Cable
Option	A pair of Red-Green glasses / A pair of Polarized glasses / RS232C / Floor Stand / Wall Bracket

Optotypes

The TSCP-800 satisfies all methods of optometry tests providing the best environment.

A type

E 0.05	9 0.05	0.05	EME 0.1 MWE 0.15	0.1
FRA 0.1 HDN 0.15	592 0.1 386 0.15	0.1 0.15	WMEW 0.2 M3WME 0.3 E3WME 0.4	0.2
KNVHO 0.2 LRF AZ 0.3 DTCEN 0.4	89259 0.2 53698 0.3 88532 0.4	0.2 0.3 0.4	EM3M 0.5 W3ME 0.6 3W3M 0.7	0.5
THFOZ 0.5 NETDP 0.6 ALVAH 0.7	26389 0.5 35926 0.6 82693 0.7	A TYPE	W3ME 0.8 EM3E 0.9 ME3M 1.0	0.8
HVNAO 0.8 LKZCE 0.9 CTFDK 1.0	88382 0.8 26953 0.9 38269 1.0	+	EM3M 1.2 W3E3 1.5 3E3M 2.0	1.2
TDVAH 1.2 LONHT 1.5 KEZFR 2.0	88635 1.2 23588 1.5 35298 2.0	+	+	+

B type

E 400	9 400	200	EME 200 MWE 150	200
FRA 200 HDN 150	592 200 386 150	100 70	WMEW 100 M3WME 80 E3WME 70	100
KNVHO 100 LRF AZ 80 DTCEN 70	89259 100 53698 80 88532 70	50 30 20	EM3M 60 W3ME 50 3W3M 40	60
THFOZ 60 NETDP 50 ALVAH 40	26389 60 35926 50 82693 40	B TYPE	W3ME 30 EM3E 25 ME3M 20	30
HVNAO 30 LKZCE 25 CTFDK 20	88382 30 26953 25 38269 20	+	EM3M 20 W3E3 15 3E3M 10	20
TDVAH 20 LONHT 15 KEZFR 10	88635 20 23588 15 35298 10	+	+	+

Specifications and design are subject to change without prior notice for improvement.



so small & yet so big



SCIENCETERA

SCIENCETERA Co., Ltd. (in Ampelos Group)
RM.404, 1st Innoplex B/D, 371-47, Gasan-dong, Geumchon-gu, Seoul, Korea 153-803
SALES DEPARTMENT (AMPELOS ENTERPRISE Co.,Ltd. in Ampelos Group)
RM.406, 1st Innoplex B/D, 371-47, Gasan-dong, Geumchon-gu, Seoul, Korea 153-803
Tel +82-2-6679-0159 Fax +82-2-6679-0147
E-mail : sales@ampelos.co.kr

DISTRIBUTED BY

SCIENCETERA