

Surgiline 900 Eye Surgery Table

Ref: M2T 400 00

**Head section positioning**

Biçaklılar Eye Surgery Table is designed for ophthalmic surgical procedures and examinations. Eye Surgery Table is constructed with high quality materials to provide hygiene and long operating life in hospitals and clinics.

Eye Surgery Table consists of head, backrest, seat and leg sections. Multifunction head section can be adjusted to required position. In order to provide larger operating field, the backrest section is in a tapering shape towards the head section. Adjustable side supports can be easily mounted to backrest section depending on the table's position of use (supine or sitting positions).

The central brake system controlled from both wheels at foot-end allows all casters to move freely, stabilize, two front casters in the straight direction, or secures the table at desired position.

Two handles on backrest section of the table provides easy maneuverability while conductivity is provided by an anti-static castor.

The tabletop can be raised and lowered, and backrest section can be moved forward and backward with remote control by hand or with a multi-functional foot pedal (optional accessory).

Three positions can be preset in memory by use of the remote control unit. For patient comfort the seat section horizontally moves concurrently with backrest section in combination with leg section movements. These positions are adjusted by an electrically operated piston.

Rechargeable battery provides uninterrupted power and can be removed easily for recharging. In urgency all the movements of table can be halted by one-touch button located on the battery control unit.

Trendelenburg position is adjusted pneumatically by a foot pedal.



Technical Details

Height adjustment: 610 - 810 mm
Table horizontal length: 1565 mm
Table base: 665 x 800 mm
Head section: 180 x 200 mm
Backrest section: 530 x 600 mm
Seat section: 530 x 540 mm
Leg section: 510 x 440 mm
Voltage: DC 24V with battery

Packing Details

Net weight: 110 kg
Gross weight: 165 kg
Packing dimensions: 1000 x 1680 x 1100 mm

Positions

Head section: Upward 100 mm
 Downward 100 mm
 Extendable 160 mm
Backrest section:
 Upward 75° yukarı
Trendelenburg: 15°
Leg section: Downward 45°



Standard Accessories

- ▶ Arm boards (2 ea)
- ▶ Flexible cover bow (1 ea)
- ▶ Remote control holder (1 ea)
- ▶ Battery recharge unit (1 ea)



Arm boards
 Ref: M9F 200 0 0



Flexible cover bow
 Ref: M2T 116 6 0



Remote control holder
 Ref: M2T 310 6 0

Battery recharge unit
 Ref: M2T 110 6 0

Accessories

- ▶ Spare battery
- ▶ Multi-functional foot pedal
- ▶ Operator chair



Spare battery
 Ref: ML5 21 001

Multi-functional foot pedal
 Ref: M2T 320 0 0

Operator Chair for Eye Surgery Table

Operator chair is designed to be used with Eye Surgery Table providing comfortable and desired seating positions. Special arm supports allow stabilization of operator's arm/hand motions. Seat height, back section and arm support positions are adjustable.



Operator chair
 Ref: M2T 150 0 0